

# VA Kidney Native iPhone App Build Screen Functionality Challenge #2 - Deployment Guide

## Revision History

Author	Revision Number	Date
TCCODER	1.0	Dec 25, 2017
TCCODER	1.1	Feb 04, 2018
TCCODER	1.2	Mar 04, 2018
TCCODER	1.3	Apr 01, 2018

## [Deployment Instructions](#)

### [1. Deployment Dependencies](#)

### [2. Organization of Submission](#)

### [3. 3rd party Libraries](#)

### [4. Configuration](#)

#### [4.1. Configuration file](#)

#### [4.2. Sample data](#)

### [5. Deployment Instructions](#)

#### [5.2. Build and run the app in a simulator or on a real device](#)

### [6. Verification](#)

### [7. Resource Contact List](#)

# Deployment Instructions

## 1. Deployment Dependencies

Before performing a deployment, it is assumed that the following have been set up:

- Xcode 9.2+
- OS X 10.12.6 or above
- iOS SDK 11 or above
- iPhone device or simulator with iOS 10+

## 2. Organization of Submission

- *src* – this directory contains the source code
- *src/VAKidneyNutrition.xcworkspace* – Xcode workspace to open.
- *docs* – this directory contains the documents for this application, including this deployment guide

## 3. 3rd party Libraries

**SwiftJSON** - <https://github.com/SwiftyJSON/SwiftyJSON>

SwiftJSON makes it easy to deal with JSON data in Swift. Version: 4.0.0

**Charts** - <https://github.com/danielgindi/Charts>

Version: 3.0.5

All libraries are configured in *src/Podfile*

## 4. Configuration

### 4.1. VAKidneyNutrition/Supporting Files/configuration.plist

You can access *configuration.plist* in Xcode in *VAKidneyNutrition.xcodeproj* project - *Supporting Files/configuration.plist*

**configuration.plist** file provides the following options:

- **ndbApiBaseUrl** - NDB base URL for API (see [USDA nutrient database](#))
- **ndbApiKey** - NDB API key (see [USDA nutrient database](#))
- **fdaApiBaseUrl** - FDA base URL for API (see [FDA Drug Interaction and Product Labeling Database](#))
- **fdaApiKey** - FDA API key (see [FDA Drug Interaction and Product Labeling Database](#))

### 4.2. Sample data

Sample data (used to fill the prototype with data) are stored in JSON files in *VAKidneyNutrition/Supporting Files/Sample Data/* group.

*allGoals.json* and *labValues.json* files added and define dependency of the generated goals and shown major lab values in Charts screen.

## 5. Deployment Instructions

### 5.2. Build and run the app in a simulator or on a real device

Pods directory should be pulled using the following command runned from `src` directory:

```
$ pod install
```

To build and run the app in a simulator or on a real device you will need to do the following:

1. Open `src/VAKidneyNutrition.xcworkspace` in Xcode
2. Select `VAKidneyNutrition` scheme from the top left drop down list.
3. Select a real iPhone (when connected) or a simulator from the top left dropdown list.
4. Click menu Product -> Run (Cmd+R)
5. Follow the verification steps in [7. Verification](#)

## 6. Verification

Follow the [challenge description](#) and [forum messages](#) to verify the app. See some notes below. Also you can follow the video (how to launch the server and verify the screens) - <https://youtu.be/CMIuq7SDJ0o> .

### Notes

- You must remove previously installed app on your device/simulator because Core Data model changed and requires fresh install.
- See details about fixed issues in [README.rtf](#)

## 7. Resource Contact List

Name	Resource Email
TCCODER	Through TopCoder Member Contact